

A Full Impression is required for each ear. It must include the full helix, crus of the helix, tragus and antitragus.

IMPORTANT: The impression must be taken just past the second bend of the ear canal.

We also request that the impression material be “compressed” at the superior junction point of the Tragus. See marked area:



- The impression should be made using a silicone base impression material. (A/B) The impression material should be medium to high viscosity material.
- The impression should be taken with an “open mouth position.”

EAR MOLD IMPRESSIONS TAKEN WITH A LIQUID POWDER MIX WILL NOT BE ACCEPTED